

YOU DECIDE THE WAY OF REPORTING AND ARCHIVING



SMALL COMPACT THERMAL PRINTER CONNECTED WITH THE FLOW METER VIA CABLE OR WIRELESSLY.

Automatic data storage onto an inbuilt SD card. SD memory card can be easily removed and inserted into your computer for subsequent storage or processing records if desired.



UROFLOWMETRY DEVICE CONNECTED TO A COMPUTER OR NOTEBOOK VIA CABLE OR WIRELESSLY.

Advanced SW functionality includes automatic/manual uroflowmetry, patient database, real time flow curves, multiple types of nomograms.

MMTD_1

MEDKONSULT medical technology s.r.o.
Pasterova 15, 779 00 Olomouc, Czech Republic

T +420 581 113 030
E mmt@mmtsystems.com



The waterproof unit allows easy cleaning and disinfection. The battery version provides incredible battery performance. No need for daily charging.



Uroflowmetry stand is easy to open/lock and adjust the height. Designed to fit underneath the commode chair.



Portable „all in case“ version with a folding stand allows easy travelling and more flexibility. It is easy to assemble and start.




Lightweight and easy to clean commode chair for kids and female patients (suitable also for videourodynamics procedures). The height adjustable MMT uroflowmetry stand is comfortable to slide underneath it.

www.mmtsystems.com



GET INTO THE FLOW WITH MMT FLOW METERS



Waterproof unit allows easy cleaning and disinfection. Battery version – no need for daily charging.


Danflow

The product is in approval process and is not placed on the market yet and is not intended to be used. Das Produkt befindet sich im Zulassungsverfahren und ist noch nicht auf dem Markt und nicht zur Verwendung vorgesehen.

MMT ACCESSORIES FOR UROFLOWMETRY MEET THE HIGH DEMANDS OF STERILE HOSPITAL ENVIRONMENT AND COMFORTABLE HANDLING DURING THE EXAMINATION



Reusable stainless steel container sterilizable in autoclave



Reusable funnel, 120°C resistant



Plastic containers with or without measuring scale



Disposable paper funnel

FLOW METERS UTILIZE WIRED / WIRELESS WEIGHT CELL TRANSDUCER TECHNOLOGY TO TRANSFORM URINE FLOW FROM THE PATIENT INTO A GRAPHICAL FORMAT



Cable connection length up to 6 meters




Unique technology at 868/915 MHz
Massively improved range between flow rate and printer unit or PC



Wireless technology at Bluetooth



WITH THE DANFLOW RANGE OF URINE FLOW METERS, MMT HAVE MADE UROFLOWMETRY SIMPLE AGAIN!



BASIC FEATURES OF FLOW METERS DANFLOW

EASY TO OPERATE
EASY TO INSTALL
METHODS CONFORM THE ICS STANDARDS

Basic features

	DANFLOW Chord CU	DANFLOW Wave CU	DANFLOW Chord PC	DANFLOW Wave
Nomograms	-	-	•	•
Real time curves	-	-	•	•
Automatic /manual mode	•	•	•	•
Patients' database	-	-	•	•
Post-processing	-	-	•	•
PDF protocol	•	•	•	•
Customized protocols	-	-	•	•
Language versions	•	•	•	•
Online remote service support	-	-	•	•
Dicom	-	-	•	•
Wired connection	•	-	•	-
Wireless connection - 868/915 MHz	-	•	-	•
Wireless connection - Bluetooth	-	•	-	•
Thermoprinter	•	•	-	-
UFM stand - standard	•	•	•	•
UFM stand - small	•	•	•	•
UFM stand - portable	•	•	•	•
Commode chair	•	•	•	•
Reusable / disposable accessories	•	•	•	•

• available - not available